

Intragam-P®								DELIVER TO			
								Weight kg			
								Dose Modified No			
								ADMINISTRATION USE			
DAY	DATE	TIME	IMMUNOGLOBULIN PRODUCT	DOSE	FREQUENCY	ROUTE	INSTRUCTIONS	DR SIGN	GIVEN BY	CHECK	TIME
			Intragam-P®	g		IV	*Infuse as per product data sheet				
			Intragam-P®	g		IV	*Infuse as per product data sheet				
			Intragam-P®	g		IV	*Infuse as per product data sheet				
			Intragam-P®	g		IV	*Infuse as per product data sheet				
			Intragam-P®	g		IV	*Infuse as per product data sheet				
			Intragam-P®	g		IV	*Infuse as per product data sheet				

Rounding: Intragam-P® is available in 3 g/50 ml and 12 g/200 ml. Round to the nearest bottle size.

Note: •Commence infusion at 1ml/minute, increasing after 15minutes gradually increase to a maximum of 3 to 4ml/minute over a further 15minutes